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CERTIFICATE

EC Certificate No. 1434-IVDD-460/2021
 EC Design-examination
 Directive 98/79/EC concerning
in vitro diagnostic medical devices

Polish Centre for Testing and Certification certifies
 that manufactured by:

Safecare Biotech (Hangzhou) Co., Ltd
 311121 Hangzhou, Building 2/203, No. 18 Haishu
 Rd.Cangqian Subdistrict, Yuhang District, China

in vitro diagnostic medical devices
 for self-testing

COVID-19 Antigen Rapid Test Device (Swab)

in terms of design documentation, comply with requirements
 of Annex III (Section 6) to Directive 98/79/EC (as amended)
 implemented into Polish law,
 as evidenced by the audit conducted by the PCBC
 Validity of the Certificate: from 23.09.2021 to 27.05.2024

The date of issue of the Certificate: 23.09.2021

The date of the first issue of the Certificate: 23.09.2021

Issued under the Contract No. MD-117/2021
 Application No: 218/2021
 Certificate bears the qualified signature.
 Warsaw, 23.09.2021
 Module A1

Vice-President

CE Zertifikat

EC Declaration of Conformity

according to the Directive 98/79/EC
 (For self-testing of ANNEX II of IVDD)

Manufacturer: Safecare Biotech (Hangzhou) Co., Ltd.
Address: Building 2/203, No.18 Haishu Rd.Cangqian Sub-district, Yuhang District, Hangzhou, Zhejiang China 311121
 Tel/Fax: +86 571 81389219 Email: admin@safecare.com.cn

EC Representative: NIC GmbH
 Erlenweg 13,49076 Osnabrück,Germany

We, the manufacturer, declare under our sole responsibility that
the medical device(s) Product Name COVID-19 Antigen Rapid Test Kit(Swab)
 For Self-testing

Type/model, identification of product allowing traceability (Where applicable) Cassette(Cov Ag-6012H)

of Category : Self-test of Annex II

is/are in conformity with the relevant provisions and requirements of Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

Applied harmonised standards, national standards or other normative documents	EN ISO23640:2015	EN ISO 18114-1:2011
	EN 13612:2002	EN ISO15223-1:2016
	EN 13641:2002	EN13532: 2002
	EN ISO 14971:2019	EN 62366:2015
	ISO13485:2016	EN ISO 17511:2003

Conformity assessment procedure **EC Declaration of Conformity(Annex III,- Section 6)**

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Signed on: 2021. 6. 10 Place: Hangzhou, Zhejiang, China

Signature (on behalf of the manufacturer) 

Name of authorized signatory:  Kebin, Qin

Position held in the company: General Manager

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Formblatt für In-vitro-Diagnostika / Form for In Vitro Diagnostic Medical Devices

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Safecare Biotech (Hangzhou) Co., Ltd

Clinical Evaluation Report

1. **Purpose:**
In order to verify the clinical performance of the improved test
2. **Material:**
Fresh negative COVID-19 samples were collected from the hospital and validated by PCR.
Fresh positive COVID-19 samples were collected from CDC and validated by PCR.
Product used: COV20082701
3. **Protocol:**
 - 3.1 **Sample Size:**
Positive Sample: >100
Negative Sample: >150
 - 3.2 **Sample's collection:**
Nasal swab specimen or nasopharyngeal swab specimen can be used by Safecare COVID-19 Antigen Rapid Test Kit(Swab) to detect the presence of SARS-CoV-2 antigen in the specimen. Internal validation studies based on Matrix Equivalency were performed on both nasal swab specimens and nasopharyngeal swab specimen, no statistic difference was observed among those specimens. All swabs were randomly blinded and assigned to testing with PCR assay as the comparator method for this study.
 - 3.3 **Sample Entry criteria:**
The samples from hospital outpatient screening cases and COVID-19 Patients who presented within 7 days of symptom onset;
Samples of people that gender and age are not limited.
 - 3.4 **Sample Exclusion criteria:**
Samples without PCR test results;
Samples that the quantity is not enough to complete the test;
Samples with failed test results (C-line has not appeared);
Freeze samples repeatedly.
 - 3.5 **Comparator method**
All samples was confirmed by PCR.
PCR tests used from Sansure Biotech Inc. and performed on ABI7500.
4. **Operator and site:**
 - Site 1:
Study Site Info: ZHEJIANG PROVINCIAL CENTER FOR DISEASE CONTROL AND PREVENTION
Researcher: Dr. ZHANG LEI
Lab Name (or Hospital or Doctor's office): Immunology Laboratory
Address: 3399 Binsheng Road, Binjiang District, Hangzhou City, Zhejiang Province
 - Site 2:
Study Site Info: THE FIRST AFFILIATED HOSPITAL ZHEJIANG UNIVERSITY SCHOOL

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COVID-19 Antigen Rapid Test Kit (Swab)

For Self-Testing Package Insert

For self testing and in vitro diagnostic use only.

INTENDED USE

The COVID-19 Antigen Rapid Test Kit (Swab) is a single-use test kit intended to detect the novel coronavirus SARS-CoV-2 that causes COVID-19. This test is intended for home use with self-collected nasal swab samples in individuals aged 18 and older. Sampling from anyone under the age of 18 should be performed under the guidance of an adult. People who are unable to carry out the test on their own should seek support. This test utilizes the lateral flow immunosassay technology for the detection of nucleocapsid protein antigen in individuals with known or suspected COVID-19. Positive results are indicative of the presence of SARS-CoV-2. Individuals who test positive should self-isolate and seek additional care from their healthcare provider. Positive results do not rule out bacterial infection or co-infection with other viruses. Negative results do not rule out SARS-CoV-2 infection. Individuals who test negative and continue to experience COVID-like symptoms should seek follow up care from their healthcare provider.

PRINCIPLE

The COVID-19 Antigen Rapid Test Kit (Swab) is an immunochromatographic membrane assay that uses highly sensitive antibodies to detect SARS-CoV-2 nucleocapsid protein from nasal swab specimens. SARS-CoV-2 specific antibodies are immobilized onto the test region of the membrane and combined with other reagents/pads to construct a test strip. During testing, the specimen reacts with anti-COVID-19 antibodies conjugated to colored particles and precoated onto the sample pad of the test. The mixture then migrates upward on the membrane chromatographically by capillary action and reacts with the reagents in test line region. Therefore, if the specimen contains COVID-19 antigen, a colored line will appear in test line. If the specimen does not contain COVID-19 antigen, no colored line will appear in the test line region, indicating a negative result. To serve as a procedural control, a colored line will always appeared in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred. The test is designed to detect nucleocapsid protein antigens in nasal swab specimens, which is different from the mutation sites have occurred in the spike protein, so it is theoretically able to detect variants including those in UK, India, South Africa and Brazil.

KIT COMPONENTS

Component	1 Test /Kit	5 Tests /Kit
COVID-19 Antigen Test	1	5
Extraction tube with buffer	1	5
Sterilized nasal swab	1	5
Package insert	1	1

ADDITIONAL SPECIAL EQUIPMENT

Timer

POTENTIAL RISKS AND BENEFITS

- Potential risks include:
 - Possible discomfort during sample collection.
 - Possible incorrect test results (see Read Results section).
 - Potential benefits include:
 - The results, along with other information, can help your healthcare provider make informed recommendations about your care.
 - The results of this test may help limit the spread of COVID-19 to your family and others in your community.

WARNINGS AND PRECAUTIONS

- Do not use after expiration date. Do not use if kit is damaged or open. Do not reuse the tests.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Discard the used testing materials in accordance with local regulations.
- The extraction buffer contains a soft sedulous agent, which in the contact contacts the skin or eyes, flush with copious amounts of water.Don't swallow the buffer.When swallowing the buffer, rinse the mouth thoroughly with water and give plenty of water to dilute the substance.If any discomfort, seek medical attention immediately.
- Children and elder please use the test under the guardian.

STORAGE AND STABILITY

Store unused test devices unopened at 4°C-30°C. If stored at 4°C-8°C, ensure that the test device is brought to room temperature before opening. The test device is stable through the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit over 30°C.

TEST PROCEDURE

Open the kit box, Check the components before use.

Please read all instructions carefully before you begin.


- Get a fat area ready, like a table. Make sure it is clear, clean and dry.
- Take off hand jewelry.
- Wash your hands for 20 seconds. Use soap and warm water, or hand sanitizer. Dry your hands using clean disposable paper towels.
- For better protection and avoid cross-contamination, disposable gloves, masks and eye protection (not provided in the package) are recommended.

[Specimen Preparation]

- Press the hole on the box, peel off the foil film on the extraction tube and insert into the hole.
- Open nasal swab package at the sticky end and take the nasal swab out.
- Gently insert the soft tip of the nasal swab into left nostril about 2.5cm (1 inch) for adult. Note: For child, the maximum depth of insertion into the nostril may be less than 2.5 cm and should be carefully and appropriately adjusted by the person, who collects sample.
- Firmly brush against the inside of the nostril in a circular motion 5 times or more.
- After the nasal swab to the right nostril and repeat the previous action. Make sure an adequate sample is collected.
- Insert the nasal swab into the tube which contains extraction buffer.
- Rotate nasal swab at least 8-10 times while pressing nasal swab tip against the bottom and the side of the tube.
- Remove the nasal swab while squeezing and turning the nasal swab against the sides of the tube to release as much liquid as possible.
- Cover the tube with dropper tip tightly and insert the tube back into the box.

[Test Procedure]

- Open the sealed pouch and take out the test cassette. For best results, the test should be performed in one hour.
- Hold the tube vertically upside down over the sample well.
- Add 3 drops specimen into the sample well by gently squeezing the sides of the tube, then start the timer.

Peel off the foil film on the extraction tube and insert into the hole

Insert nasal swab to an appropriate depth, not less than 2.5 cm and rotate 5 times in each nostril

Rotate the swab 8-10 times

Press the swab and release as much liquid as possible

Place the cap

Add 3 drops to sample well(s)

4.Wait for colored lines to appear.The test result can be read in 10-15 minutes. DO NOT read after 20 minutes.

[After the testing]

1. After you have done the test, put all parts of the kit in the waste bag. Discard the waste bag in accordance with local regulations.

2. If you are doing more than 1 test, clean the table with 75% alcohol or sanitizer. Wash your hands between each test.

INTERPRETATION OF TEST RESULT

POSITIVE RESULT		A colored line appears in the control line region (C) and a colored line appears in test line region (T). <p>*NOTE: The intensity of the color in the test line region will vary depending on the concentration of COVID-19 antigen in the specimen. Therefore, any shade of color in the test line region should be considered positive.</p>
NEGATIVE RESULT		A colored line appears in the control line region (C) and no line appears in test line region (T).
INVALID RESULT		No line appears in control line region(C). Insufficient buffer volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the procedure with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

BUILT-IN CONTROL

This test contains a built-in control feature, the C line on the test. The C line develops after adding sample solution. Otherwise, review the whole procedure and repeat test with a new device.

WHAT SHOULD I DO AFTER TEST

If the test result is positive

- There is currently a suspicion of a COVID-19 infection
 - Contact your doctor / general practitioner or the local health department immediately
 - Comply with local guidelines for self-isolation
 - To have a PCR confirmatory test performed

If the test result is negative

- Continue to comply with all applicable rules regarding contact with others and protective measures
 - An infection may also be present if the test is negative
 - If it is suspected, repeat the test after 1 - 2 days, as the coronavirus cannot be accurately detected in all phases of an infection

If the test result is invalid

- Possibly caused by incorrect test execution
 - Repeat the test
 - If the test results remain invalid, contact a doctor or a COVID-19 test center

Note: Do not take any decision of medical relevance without first consulting your medical practitioner

PERFORMANCE CHARACTERISTICS

- Limit of Detection: The Limit of Detection (LoD) was determined by evaluating different dilutions of heat inactivated SARS-CoV-2 virus. LoD is defined as the virus concentration - 1.3 x 10⁶ TCID₅₀/mL at which a minimum of 19 replicates out of 20 generate a Reactive result.
- High-dose Hook Effect: No hook-effect was observed at 3.33X10⁷ TCID₅₀/mL of inactivated SARS-CoV-2 virus culture.
- Clinical study: A comparison was conducted with RT-PCR, according the test data listed below table, the Relative Sensitivity is 96.4% (212/220), the Relative Specificity is 100%(200/200) and the overall agreement is 98.1% ((212+200)/420).

		PCR result		Total
		Positive	Negative	
Safecare Test	Positive	212	0	212
	Negative	8	200	208
	Total	220	200	420

4. Cross-reactivity: Cross-reactivity studies were performed to demonstrate that the test does not react with the following microorganisms in the table below at concentration of 1x10⁷ TCID₅₀/mL for viruses and 1x10⁸ CFU/mL for bacteria.

Human metapneumovirus (hMPV)	Human parainfluenza virus 1	Adenovirus	<i>Bordetella pertussis</i>	<i>Streptococcus pneumoniae</i>
Human coronavirus OC43	Human parainfluenza virus 2	Rhinovirus	<i>Chlamydia pneumoniae</i>	<i>Streptococcus pyogenes</i>
Human coronavirus 229E	Human parainfluenza virus 3	Enterovirus	<i>Haemophilus influenzae</i>	<i>Mycobacterium tuberculosis</i>
Human coronavirus NL63	Human parainfluenza virus 4	Influenza A	<i>Legionella pneumophila</i>	<i>Staphylococcus aureus</i>
Respiratory Syncytial Virus	MERS	Influenza B	<i>Mycoplasma pneumoniae</i>	<i>Candida albicans</i>

5. Interference: The following endogenous interference substances were evaluated at the concentrations listed and no effect was found.

Whole blood (2%), three OTC nasal sprays (10%), three OTC nasal drop (25%), three nasal mouthwashes (25%), 4-Acetamidophenol (10mg/mL), Acetylsalicylic acid (20mg/mL), Chlorpheniramine (5 mg/mL), Dextromethorphan (10mg/mL), Diphenhydramine (5mg/mL), Ephedrine (20mg/mL), Guaifolins glycol ether (20mg/mL), Oxymetazolin (100mg/mL), Phenylephrine (100mg/mL), Phenylephrandiamine (20mg/mL), Oseltamivir Phosphate (10mg/mL), Mupirocin (10mg/mL), Vitamin A (10%), D-Pantened (10%)

LIMITATIONS AND POSSIBLE ERRORS

- The COVID-19 Antigen Rapid Test Kit (Swab) is intended for use as a self-test and may only be used for qualitative detection of SARS-CoV-2 antigens. The colour intensity of a positive line shall not be evaluated as quantitative or semi-quantitative.
- The COVID-19 Antigen Rapid Test Kit (Swab) should be only for the detection of SARS-CoV-2 antigen, not for any other viruses or pathogens.
- The performance was evaluated using the procedures provided in this product insert only. Modifications to these procedures may alter the performance of the test.
- A negative result does not exclude the possibility of COVID-19 infection.
- The results obtained with this assay, especially in the case of weak test lines which are difficult to interpret, should retest or go to a medical institution for testing.
- The test is intended for infection detection and not for determine infection status.The test is used for the auxiliary diagnosis of patients with the COVID-19 and can not be used as the sole diagnostic indicator of whether the test subject is infected with the COVID-19.

ACCESSORIES

Accessory	Manufacturer	EC Representative	CE Mark
Swab A	Jiangsu Changfeng Medical Industry Co.,Ltd. <p>Toujiao Town,Guangding District <p>Yangzhou Jiangsu 225109, P.R.China</p></p>	Lins Service & Consulting GmbH <p>Obere Seegasse 34/2,69124 Heidelberg <p>Germany</p></p>	
Swab B	Goodwood Medical Care Ltd. <p>1-2 Floor,3-919 Yongzheng Street, Jinzhou District,Dalian 116100 Liaoning, P.R.China</p>	CMC Medical Devices & Drugs S.L. <p>CHoracio Lengo <p>No.18,CP29006,Málaga,Spain</p></p>	
Swab C	QIOTEST LABWARE MANUFACTURING CO.,Ltd. <p>No.339 Beihai West Road,Haimen, 226100 Jiangsu, P.R. China</p>	Wellking Ltd <p>Enterprise Hub,NW Business Complex,1 Beraghmore Road,Derry,BT48 8SE <p>Northern Ireland</p></p>	
Swab D	Jiangsu HanHeng Medical Technology Co.,Ltd. <p>16-B4,#1 North Qingyang Road, Tiaanning District, Changzhou, 213017 Jiangsu P.R. China</p>	Luxus lebenswelt GmbH <p>Koshstr.1,47877,Willich,Germany</p>	
Swab E	Jiangsu Rongye Technology Co.,Ltd <p>Toujiao Town,Yangzhou City,225109 Jiangsu, P.R.China</p>	Riomavix S.L. <p>Calle de Almirans 55,1D, Madrid 28039 Spain</p>	

INDEX OF SYMBOLS

	Do not reuse		In vitro diagnostic use only
	Stored between 4-30°C		Consult instruction for use
	Keep away from sunlight		Lot number
	Use-by date		Contains sufficient for <-> tests
	Manufacturer		Date of Manufacture
	Sterilized using ethylene oxide		Sterilized using irradiation
	Authorized Representative in the European Community		CE Mark
	Importer		Safecare Biotech(Hangzhou) Co., Ltd. <p>Building 2/203, No.18 Haishu Rd., Cangjiang Sub-district Yuhang District, Hangzhou, 311121, China <p>Tel/Fax: +86 571 81389219 <p>Email: admin@safecare.com.cn</p></p></p>
	NIC GmbH <p>Erlenweg 13, 49076 Osnabrück, Germany <p>Mobil: +49 (0) 9116706 <p>Email:info@nic-industry.com</p></p></p>		EC REP

Angela G. GmbH
Bahnhofstr. 2-6, 66953 Pirmasens DE
Tel.: +49 (0) 69 95635730
Mobil: +49 (0) 176 3117 7329
Email: info@angela.de

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Accessorio	Produttore	Rappresentante CE	Marchio CE
Tampone A	Jiangsu Changfeng Medical Industry Co.,Ltd. Toujiao Town,Guangling District Yangzhou Jiangsu 225109, P.R.China	Ulins Service & Consulting GmbH Obere Seegasse 34/2,69124 Heidelberg Germany	
Tampone B	Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yongzheng Street,Jinzhou District,Dalian 116100 Liaoning, P.R.China	CMC Medical Devices & Drugs S.L. C/Horacio Lengo No.18, CP29006,Málaga, Spain	
Tampone C	CITOTEST LABWARE MANUFACTURING CO.,Ltd. No.339 Beihai West Road,Haimen, 226100 Jiangsu, P.R. China	Wellkang Ltd Enterprise Hub,NW Business Complex,1 Berghmore Road,Derry,BT48 8SE Northern Ireland	
Tampone D	Jiangsu Han-Heng Medical Technology Co.Ltd. 16-B4-#1 North Qingyang Road, Tianning District, Changzhou, 213017 Jiangsu P.R. China	Luxus lebenswelt GmbH Koshstr.1,47877,Willich,Germany	
Tampone E	Jiangsu Rongye Technology Co.,Ltd Toujiao Town,Yangzhou City,225109 Jiangsu, P.R.China	Riomavix S.L. Calle de Almansa 55, 1D, Madrid 28039 Spain	

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Date: 29 Nov. 2021

To whom it may concern,

We, **Safecare Biothech (Hangzhou) Co., Ltd.**, having our office at F1.2 Blog.2, No.18 Haishu Road, Hangzhou 311121 China, as a manufacturer of COVID-19 Antigen Rapid Test Kit(Swab) and COVID-19 Antigen Rapid Test Kit(Saliva), hereby declare that our COVID-19 Antigen Rapid Test Kit(Swab) and COVID-19 Antigen Rapid Test Kit(Saliva) remain effective for the detection of SARS-CoV-2 antigen even in the emergence of newly discovered variant **Omicron** found in South Africa.

Since the recognition site of the raw materials used in our antigen test is the nucleocapsid protein (nucleoprotein or protein N) antigens, which is different from the mutation sites, we would like to declare that our products are theoretically able to detect this new variant in South Africa.

Lastly, We Safecare will strictly implement our quality management system and strive to provide the best products to the customers. We will also inform you officially if there is any update information of our COVID-19 Antigen Rapid Test Kit(Swab) and COVID-19 Antigen Rapid Test Kit(Saliva).

Yours sincerely,

SAFECARE BIOTECH (HANGZHOU) CO., LTD.
HANGZHOU SAFECARE BIOTECH CO.,LTD.

EC Declaration of Conformity

according to the Directive 98/79/EC
(For self-testing of ANNEX II of IVDD)

Manufacturer: Safecare Biotech (Hangzhou) Co., Ltd.

Address: Building 2/203, No.18 Haishu Rd.Cangqian Sub-district, Yuhang District, Hangzhou, Zhejiang China 311121
Tel/Fax: +86 571 81389219 Email: admin@safecare.com.cn

EC Representative: NIC GmbH
Erlenweg 13,49076 Osnabrück,Germany

We, the manufacturer, declare under our sole responsibility that

the medical device(s)	Product Name	COVID-19 Antigen Rapid Test Kit(Swab) For Self-testing
	Type/model, identification of product allowing traceability (Where applicable)	Cassette(Cov Ag-6012H)
of Category	: Self-test of Annex II	

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Conformity assessment procedure EC Declaration of Conformity(Annex III,- Section 6)

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Signed on: 2021.6.10 **Place:** Hangzhou, Zhejiang, China

Signature (on behalf of the manufacturer) _____

Name of authorized signatory: Kebin, Qiu
Position held in the company: General Manager
Seal/Stamp:

CERTIFICATE

EC Certificate No. 1434-IVDD-460/2021

**EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices**

Polish Centre for Testing and Certification certifies
that manufactured by:

**Safecare Biotech (Hangzhou) Co., Ltd
311121 Hangzhou, Building 2/203, No. 18 Haishu
Rd.Cangqian Subdistrict, Yuhang District, China**

in vitro diagnostic medical devices
for self-testing

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in terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC
Validity of the Certificate: from 23.09.2021 to 27.05.2024

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Issued under the Contract No. MD-117/2021
Application No: 218/2021
Certificate bears the qualified signature.
Warsaw, 23.09.2021
Module A1

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

Safecare Covid-19 Antigen Rapid Test (Swab)

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COVID-19 In Vitro Diagnostic Devices and Test Methods Database

COVID-19 In Vitro Diagnostic Medical Device - detail

COVID-19 Antigen Rapid Test Kit (Swab)

Manufactured by Safecare Biotech (Hangzhou) Co. Ltd, China - www.safecare.com.cn

Device identification number	1489
CE Marking	Yes
HSC common list	Yes
HSC mutual recognition	Yes
Format	Near POC / POC
Target	Antigen
Specimen	Nasal swab
Commercial Status	Commercialised
Last Update	2021-05-12 12:58:25 CET

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Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization

**Safecare Biotech (Hangzhou)
Co., Ltd.**
Building 2/203, No. 18 Haishu Rd.
Cangqian Sub-district, Yuhang District
Hangzhou
311121 Zhejiang
P.R. China

has established and applies a quality management system for medical devices
for the following scope:

**Design and Development, Manufacture and
Distribution of In Vitro Diagnosis of
Rapid Test of Fertility, Drug of Abuse,
Cardiac Markers, Infectious Diseases**

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2020-08-02
Certificate Registration No.: SX 60149068 0001
An audit was performed. Report No.: 15096152 005
This Certificate is valid until: 2023-06-06

Certification Body

Date 2020-08-02

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-mail: cert-validity@de.tuv.com http://www.tuv.com/safety

CE

CE

EC Declaration of Conformity

according to the Directive 98/79/EC

(applicable to IVD Devices of NOT Annex II and NOT self-test)

Manufacturer: Safecare Biotech (Hangzhou) Co., Ltd.

Address: Building 2/203, No.18 Haishu Rd.Cangqian Sub-district,
Yuhang District, Hangzhou, Zhejiang China 311121

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the medical device(s) Product Name COVID-19 Antigen Rapid Test Kit(Swab)

Type/model, identification of product allowing traceability (Where applicable) Cassette(COV Ag-6012)

of Category: Common/Others IVD
(Devices of NOT Annex II and NOT self-test)

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Conformity assessment procedure Module A (EC Declaration of Conformity) (Annex III, except point 6)

Notified Body (name & number) **NOT applicable**
Certificate & number

Signed on 28th Sep.,2020 Place: Hangzhou, Zhejiang, China

Signature (on behalf of the manufacturer) Kevin Qiu 2020.9.28

Name of authorized signatory: Kevin Qiu

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- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. Trends Microbiol 2016;24:490-502.
- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019;17:181-192.
- Wong G, Liu W, Liu Y, Zhou B, Bi Y, Gao GF, MERS, SARS, and Ebola: the role of super-spreaders in infectious disease. Cell Host Microbe 2015;18:398-401.

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Safecare Biotech (Hangzhou) Co., Ltd.
Building 2/203, No.18 Haishu Rd, Cangqian Sub-district
Yuhang District Hangzhou, 311121, China

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